




Algocyte® - Quick guide

3 simple steps to get you started

Step 1 – Registration and logging in

Visit <https://uk.algocyte.com> (U.K.) or <https://eu.algocyte.com> (E.U.), to sign up, create and verify your account*.




*Refer to Step 5 of the “User Sign-up and Registration” section in the user manual to learn more about how to request your provider account access.

Step 2 – Getting Started

Once logged in, confirm you have read the T&C's before closing and proceed to complete your account


Click on **'Profile'** and ensure all your **'Demographic information'** is filled out completely



*Providers will be taken to the 'My Subscribers' window, where you will need to select a pre-linked patient.

Now it is time to order your first home kit

Click on **'Kit Management'**, you will be presented with two options: **Register Kit, Order kit**



Click on **'Order Kit'** and fill out the following

Order Kit	Register Kit
To schedule a test kit dispatch, you will need to choose: 1. The type of test kit 2. The dispatch frequency	Once your kit has arrived you should register it using either: 1. QR Code (Mobile only) 2. Enter the number manually
Test kit dispatch scheduling New Dispatch Request Select test Full blood count CRP Calcium Urea and Electrolytes LFT (not AST) Anaemia Profile - FBC, ESR, Fe, TIBC, Ferritin, B12 (active) and Folate (RBC) Step 1 Select frequency of kit dispatch ASAP Single test Daily Weekly Fortnightly Every 4 weeks Step 2	Please scan or upload your Kit order code to register it Please scan QR code or insert order number manually Order number Upload to Server Scan QR Code CLICK TO SCAN



Step 3 – My Data

Now you registered your kit and sent back your sample, its time to view your data!

The Dashboard

This is where you can access and view your health data and personalised analysis in different forms

Algocyte® Score is a numerical summary of the latest blood test results presented in a single view

Covid Profile represents the COVID-19 antibody test result, scale from 0 to 10

Personalised Test Score

Click on the 'Dashboard' tab

Results panel displays the most recent blood test result, as a bar-chart and reference range for each analyte.

'Adaptive On' - uses historical test results as a baseline to provide a more personalised analysis.

Icon

Action

- Bar chart icon: Navigates to the Trends Graphs
- Menu icon: Navigates to the Tabular Report
- Line graph icon: Navigates to the Individual Report
- Folding plot icon: Navigates to the Folding Plot
- PDF icon: Generates an Algocyte® Report
- Adaptive toggle: Turns on/off the Adaptive feature. *Note: the Adaptive feature is on by default.*

Algocyte Ring


Hovering over the ring and its different colours (Green, Amber and Red) opens more detailed tables on the data presented like the examples below.



Amber Section



Score Example



Green Section

You are now all set up for Algocyte®, the AI that learns from and monitors your health.

Please refer to the user manual for full user experience and support